

SALTER LABS®



Salter-Style®

Micro Cannula

A Salter-Style, adult, low flow cannula for use with up to 3 LPM of oxygen. Provides superior comfort for the active, low flow user.



Features

- ☐ Unique patented, anatomically correct design
- ☐ Curved, tapered nasal prongs
- ☐ Three-channel, crush-resistant safety tubing
- ☐ The “diminutive” streamlined design and clear material is almost invisible on the skin

Benefits

- ☐ Provides a natural fit for better comfort
- ☐ Better anatomical fit for increased comfort
- ☐ Uninterrupted flow of oxygen
- ☐ Feather-light weight provides comfort and reduces behind-the-ear pressure




slc-01

Salter-Style® Micro Cannula

Patient compliance can be difficult with low flow oxygen therapy users often because they are still in the mainstream of life and are self conscious about wearing a cannula. The diminutive size of the Micro Cannula makes it less visible on the face which can increase patient confidence and compliance. In addition, the Micro Cannula is soft and comfortable to wear because it is made with the *Salter-Style* manufacturing process.

Clear, lightweight cannula for oxygen flows up to 3 LPM. Unique one piece, anatomical design provides maximum patient comfort by eliminating the flap which can cause irritation. Over-the-ear style head set provides secure positioning, allowing freedom of movement.

*Curved,
tapered
nasal
prongs*



Ordering Information

Micro Cannula

Salter-Style Micro Cannula

Micro cannula with 7' supply tube (for up to 3LPM of oxygen)

50

1616

Micro cannula with 7' supply tube and female thread grip connector

50

1616-TG

Adult Micro cannula with 1' supply tube

50

1616-1

Adult Micro cannula with 4' supply tube

50

1616-4



100 W. Sycamore Road - Arvin, CA 93203

Phone: 661-854-3166

1-800-421-0024

1-800-235-4203

Fax: 661-854-3850

Toll Free Fax: 1-800-628-4690

Fax: 661-854-3850

E-mail: salterlabs@us.salterlabs.com

www.salterlabs.com

CE

0482

ISO 9001 & ISO 13485