



## Pulmo-Aide® Compact Compressor Nebulizer System

### Quality, Reliability, and Affordability

The DeVilbiss Pulmo-Aide Compact Compressor Nebulizer System combines **effective treatment with a small foot print design**. Built to last, the Pulmo-Aide Compact has metal mechanical parts, making it durable, and a **five-year warranty** ensuring quality therapy for an extended period of time.

For patients on the go, the Pulmo-Aide Compact has a convenient carrying handle and a **small, compact design**, and can be **easily transported for AC use**. Reusable nebulizers can be included for convenience, and reliable pressure and flow make it one of the highest quality units on the market.



Pulmo-Aide Compact Compressor Nebulizer System (Model 3655D)

### Performance and Design Features

- Reliable pressure and flow
- 5-year warranty
- Reusable nebulizer can be included
- Small footprint
- Small, compact design
- Convenient carrying handle
- Metal mechanical parts

## Parts and Accessories

DEVILBISS PART #	PRODUCT DESCRIPTION	DEVILBISS PART #	PRODUCT DESCRIPTION
<b>Aerosol Masks</b>		<b>VixOne Disposable Nebulizer with PediNeb™</b>	
AAM12	Adult Aerosol Mask (Package of 12)	3655PN-612	Case of 12
AAM20	Adult Aerosol Mask (Package of 50)	<b>Peak Flow Meters</b>	
DL1050	Character Aerosol Mask (Package of 50)	600-112	Case of 12
PAM10	Pediatric Aerosol Mask (Package of 50)	<b>Air Inlet Filters</b>	
<b>Disposable Nebulizers</b>		3655D-601	Package of 5
3655D-621	VixOne™ Disposable Nebulizer (Case of 50)	3655D-602	Package of 25
<b>Reusable Nebulizers</b>		<b>Carrying Case for Compact Compressor Nebulizer</b>	
900D-610	VixOne (Case of 10)	3655D-603	Each

### Technical Specifications

Built to meet the requirements of HCPCS code:	E0570
Product Weight:	3.8 lbs.
Product Width:	7.5 in.
Product Length/Depth:	7.2 in.
Product Height:	4.0 in.
Shipping Dimensions (box 1):	1 pk. 9.5"w/ 9"d/ 5"h
Shipping Dimensions (box 2):	6 pk. 19.75"w/ 15"d/ 10"h
Ship Method(s):	Ground, Next Day, 2 Day or 3 Day, Freight
Shipping Weight:	1 pk. (5 lbs.) 6 pk. (33 lbs.)
Electrical requirements:	2.5 amps, 115 VAC 60 Hz
Sound level*:	58-62 dBA
Power consumption:	140 watts maximum
Maximum compressor pressure*:	35 psig
Compressor free air flow*:	13 lpm
Operating temperature range:	40 to 104 deg. F (5 to 40 deg. C)
Supplied nebulizer operating pressure*:	12 psig (min)
Supplied nebulizer liter flow*:	5.5 lpm (min)



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA.  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)