

SALTER LABS®



Salter-Style®

Small Cannulas

The Small Salter-Style Cannulas are for oxygen flows up to 3 LPM or 6 LPM. A full range of sizes for Premature, Neonate, Infant, Intermediate Infant, Pediatric and Small Adult patients. The clear, lightweight cannula is made with the Salter-Style manufacturing process, producing a soft, anatomically designed cannula that is gentle to delicate facial skin.



Features

- ☐ Anatomically curved to fit upper lip
- ☐ Curved, tapered nasal prongs
- ☐ Over-the ear-style head set
- ☐ Three-channel safety tubing
- ☐ Soft, lightweight, clear material

Benefits

- ☐ Eliminates flap or ridge which can cause irritation
- ☐ Improved anatomical fit for increased comfort
- ☐ Provides secure positioning
- ☐ Allows uninterrupted flow of oxygen
- ☐ Helps to provide maximum comfort with less visibility on the face




slc-04

Salter-Style® Small Cannulas: Premature, Neonate, Infant, Intermediate Infant, Pediatric and Small Adult Sizes

1610-Premature... less than term
 1611-Neonate...more than 1400 grams (approx.)
 1601-Infant... newborn, full term
 1615-Intermediate Infant... 3 to 12 months (approx.)
 1602-Pediatric... 1 year or older (approx.)
 1602-Small Adult

Clear, lightweight cannula with unique, one piece anatomical design, eliminates the need for irritating flap, resulting in maximum patient comfort. Over-the-ear style headset provides secure positioning, allowing freedom of movement. Provided with three channel safety tubing.

*Soft, clear material—
anatomically
curved for
improved
upper lip fit.*



Ordering Information

*Salter-Style Small Cannulas**

	Units (per case)	Part Number		Maximum LPM
		Safety Channel	Smooth Bore	
Premature clear cannula with 7' supply tube	50	1610		3 LPM
Half case packaging of 1610 Premature cannula	25	1610HC		
Neonate clear cannula with 7' supply tube	50	1611		3 LPM
Half case packaging of 1611 Neonate cannula	25	1611HC		
Infant clear cannula with 7' supply tube	50	1601		3 LPM
Half case packaging of 1601 Infant cannula	25	1601HC		
Infant Quiet clear cannula with 7' supply tube **	50	1601Q	1651	6 LPM
Intermediate Infant clear cannula with 7' supply tube	50	1615		3 LPM
Half case packaging of 1615 Infant cannula	25	1615HC		
Pediatric, Small Adult clear cannula with 7' supply	50	1602	1652	6 LPM
Pediatric, Small Adult clear cannula with 7' supply tube and female thread grip connector	50	1602TG		6 LPM
Pediatric, Small Adult clear TLCannula™ with foam ear cushions and 7' supply tube	25	1602TLC		6 LPM

***Special lengths available*



100 W. Sycamore Road - Arvin, CA 93203
 Phone: 661-854-3166
 1-800-421-0024
 1-800-235-4203
 Fax: 661-854-3850
 Toll Free Fax: 1-800-628-4690
 Fax: 661-854-3850
 E-mail: salterlabs@us.salterlabs.com
 www.salterlabs.com

CE
0482

ISO 9001 & ISO 13485