


SALTER LABS®



Salter Labs

Tubing Connectors

Easy to use connectors for coupling oxygen tubing. Durable, long lasting and reliable.



Features

Tubing Connector 1215

- ☐ Rounded and tapered edges
- ☐ Ribbed "Sure Grip" body
- ☐ Unique multiple concentric rings
- ☐ Universal conical tips
- ☐ Vinyl "See Through" material

Benefits

- ☐ Prevents snagging or catching on various surfaces
- ☐ Easy to grasp or hold when making connections
- ☐ Ensures secure leak free coupling
- ☐ Makes connecting easy
- ☐ Can visually check for patency. Long lasting durable construction

Tubing Connector 2005

- ☐ Lightly tinted vinyl material
- ☐ Durable, long lasting. Can visually check for occlusions


Variety of uses. . . Hospital, Skilled Nursing Facility, Homecare, Clinic




slc-13

Tubing Connector 1215*

Unique oxygen tubing connector constructed of vinyl material that allows visual inspection for patency. "Ribbed Body" design allows sure grip. Rounded and tapered shoulders help prevent snagging or catching on various surfaces. Special Concentric Rings and Conical Tip design make connection easy and positive.



"COUPLED" TOGETHER



Tubing Connector 2005

A simple tubing connector.
Durable, reliable to use.
Cost effective



Ordering Information

Salter Tubing Connectors

	Units (per case)	Part Number
Oxygen tubing connector with concentric rings for easy "twist & pull" connection for connecting additional lengths of oxygen tubing.	50	1215
2" Oxygen tubing connector - For connecting additional lengths of oxygen tubing (Packages of 50 each).	N/A	2005



100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
Fax: 661-854-3850
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

CE
0482

*U.S. Patent No. 5,573,280