



### Every Unit Comes With:

- Single battery
- Nasal cannula
- AC charger for home
- DC charger for car
- Carry bag with padded adjustable strap

### Optional Accessories:

- Double Battery:  
Up to 12.75 hours of power  
Model No. 1400-3010-16
- Desktop Battery Charger  
Model No. 1400-3030
- Backpack  
Model No. 1170-3420
- Custom Cart  
Model No. 1170-3442

Getting your  
**OxyGo NEXT<sup>®</sup>** is easy.

Simply Contact Your Provider:

Introducing the  
**NEXT** generation



The new OxyGo NEXT<sup>®</sup>  
keeps delivering —  
so you can **Keep Going.**<sup>sm</sup>

Connect With Us



[oxygo.life](http://oxygo.life)

PN B156-0226

[oxygo.life](http://oxygo.life)

# Keep Going.<sup>sm</sup>



**12+**  
Hour  
Battery!\*

Breathe Easy. Live Fully.

**OxyGO NEXT**®

The **NEW** OxyGo NEXT®  
is designed to keep up with  
your life – not slow it down.

With six intelligent pulse dose settings and extended battery options, the new OxyGo NEXT delivers powerful oxygen therapy in a compact, lightweight design. Quiet, reliable, and easy to use, it gives you the freedom to keep moving, stay connected, and **Keep Going**.

OxyGo® portable oxygen concentrators (POCs) are lightweight medical devices that make doctor-prescribed medical oxygen from the air around you. You never need to worry about running out of oxygen or carrying heavy tanks again.

- Simple controls, bright digital display
- User-replaceable batteries & components
- Designed for 24/7 use, at home or away
- Meets FAA-acceptance criteria for use onboard aircraft
- Optional Bluetooth connectivity (U.S. and Canada)

Model No. 1400-3000NG

LEARN MORE: [oxygo.life](http://oxygo.life)

<b>Size, Weight</b>	<p><b>With Single Battery:</b> (L) 7.18" x (W) 3.27" x (H) 8.14" 4.8 lbs. / 2.2 kg.</p> <p><b>With Double Battery:</b> (L) 7.18" x (W) 3.27" x (H) 9.02" 5.8 lbs. / 2.6 kg.</p>
<b>Oxygen Flow</b>	<p><b>Pulse Dose Settings:</b></p> <ol style="list-style-type: none"> <li>1 (210ml/min)</li> <li>2 (420ml/min)</li> <li>3 (630ml/min)</li> <li>4 (840ml/min)</li> <li>5 (1,050ml/min)</li> <li>6 (1,260ml/min)</li> </ol>
<b>Power</b>	<p><b>AC Power Supply:</b> 100-240 V, 50-60 Hz (auto-sensing for use worldwide)</p> <p><b>DC Power:</b> Car / Boat / RV compatible</p>
<b>Recharge Times</b>	<p><b>Single Battery:</b> Up to 3 hours</p> <p><b>Double Battery:</b> Up to 4 hours</p>
<b>Battery Life</b>	<p><b>Single Battery:</b> Up to 6.5 hours</p> <p><b>Double Battery:</b> Up to 12.75 hours*</p>
<b>Noise</b>	<p>39 dBa (at setting 2) Quieter than a library!</p>
<b>Warranty</b>	<p>5-year warranty See OxyGo Warranty Statement for details.</p>

\* Double battery at flow setting 1